## Determine volumetrically, the value of 'n' i.e., the number of water molecules of crystallisation in washing soda *Na*<sub>2</sub>*CO*<sub>3</sub>.*nH*<sub>2</sub>*O*, 7.0 g of which have been dissolved per litre of the given solution. Provided M/10 HCl

**Chemical Equation:** 

 $\mathrm{Na_2CO_3.}n\mathrm{H_2O}(aq) + 2\mathrm{HCl}(aq) \longrightarrow 2\mathrm{NaCl}(aq) + (n+1)\mathrm{H_2O}(l) + \mathrm{CO_2}(g)$ 

**Indicator.** Methyl orange. **End Point.** Yellow to light pink (Acid in burette)

## **Procedure**

- 1. Rinse and fill the burette with the standard HCl solution.
- 2. Rinse the pipette with the sodium carbonate solution and pipette out 20 ml of this solution in the titration flask.
- 3. Add 2-3 drops of methyl orange indicator to the titration flask. The colour of the solution becomes yellow.
- 4. Note the initial reading of the burette and run acid solution slowly in the titration flask till the solution becomes light pink.
- 5. Note the final reading of the burette and find out the volume of acid used.
- 6. Repeat the procedure 4-5 times to get a set of at least three concordant readings.

## **Observations**

Molarity of given HCl solution =  $\frac{M}{10}$ 

Volume of sodium carbonate solution taken for each titration = 20.0 ml

S.No.	Initial reading of the burette	Final reading of the burette	Volume of the acid solution used
1.	_	-	— ml
2.	-	_	— ml
3.	-	_	— ml
4.	_	_	— ml

Concordant volume = x ml (say).

Calculations

The molarity of standard HCl = 0.1 M

Since in the balanced equation two moles of HC1 react with one mole of  $Na_2C0_3.nH_20$ .

 $\frac{M_{H\,Cl}\times V_{HCl}}{M_{Na_2CO_3}\times V_{Na_2CO_3}} == \frac{2}{1}$ 

 $\frac{0.1\times x}{M_{Na_2CO_3}\times 20.0}=\frac{2}{1}$ 

$$M_{Na_2CO_3} = \frac{1 \times x \times 0.1}{2 \times 20.0} = \frac{x}{400}$$

 $M_{Na_2CO_3} = \frac{Streanth \ per \ litre}{Molar \ mass \ of \ Na_2CO_3.nH_2O}$ 

Molar mass of  $Na_2CO_3.nH_2O = \frac{Streamth per litre}{M_{Na_2CO_3}}$ 

$$=\frac{7.0}{\frac{x}{400}}$$
 gmol<sup>-1</sup>

But molar mass of hydrated sodium carbonate = (106 + 18n) g mol<sup>-1</sup>

therefore,

 $\frac{7.0}{\frac{x}{400}} = 108 + 18n$ 

Knowing the titre value, x, the value of n can be calculated.

## Result

The number of the molecules of water of crystallisation in washing soda is ..... Note. The result is expressed as nearest whole number. The value of n is 10 in this case.