To Prepare a Pure Sample Of the Complex Potassium Trioxalatoferrate (III), Kg[Fe(C₂O₄)₃I . 3H₂0

Theory

The complex potassium Trioxalatoferrate (III) can be prepared by dissolving freshly prepared ferric hydroxide in a solution of potassium oxalate and oxalic acid.

Requirements

Three beakers (250 mL), china dish, funnel, funnel-stand, glass-rod, wash bottle, tripod stand and wire gauze. Ferric chloride, oxalic acid hydrated, potassium oxalate and potassium hydroxide.

Procedure

- 1. Dissolve 3.5 g of anhydrous ferric chloride 50 mL of distilled water in a 250 mL beaker.
- 2. In another beaker dissolve 4 g of potassium hydroxide in 50 mL of water.
- 3. Add KOH solution to FeCl₃ solution in small portions with constant stirring. Filter the precipitates of ferric hydroxide so formed through a buchner funnel. Wash the ppt. with distilled water.
- 4. In another beaker (250 mL) take 4 g of hydrated oxalic acid and 5.5 g of hydrated potassium oxalate. Add about 100 mL of water and stir thoroughly to get a clear solution.
- 5. Add the freshly prepared Fe(OH)₃ ppt. in small amounts to the above solution with constant stirring. The ppt. get dissolved. If ppt. does not dissolve then warm it and leave the contents for some time.
- 6. Filter and transfer the filtrate to china dish and heat on a sand bath or wire-gauze to obtain crystallisation point.

- 7. Now place the china dish on a beaker full of cold water and keep it aside for crystallisation. China dish should be covered with a black paper as the complex is sensitive to light.
- 8. Decant off the mother liquor, wash the crystals with a small amount of ethyl alcohol and dry them between the folds of filter paper.
- 9. Find out the weight of the crystals.

Observations

Weight of the crystals obtained =g
Colour of the crystals is......

Precautions

- 1. Do not concentrate the solution too much.
- 2. Let the concentrated solution cool slowly and undisturbed to get large crystals.